



Firm mattress consisting of a lamination of high-resilience foam and soft foam ensuring a soft, uniform feeling of comfort. Easy maintenance, the Series V-1100 allows a variety of formats. Flexible polyurethane-polyester fabric helps to reduce pressure.



Supports a weight up to 318 kg



Easy ajustable with removable side booster 91 cm to 106 cm to 121 cm



Antibacterial, antimold and antifungal

Features

- › The V-1100 Series is covered with an Integrity 30 fire net to make it self-extinguishing.
- › Meets flammability standards CAN/ULC-S137:2017, CAL TB-129 and CAL 117.
- › Removable booster allow you to adjust the size of the mattress if necessary.
- › Zippered cover made of antibacterial, antifungal, antimold and fire resistant polyurethane laminated polyester.
- › 15 cm height.



VISCO-ELASTIC FOAM
High-density 64 kg/m³
5 KGF compressive

High Density HR Foam
Ultra Compressed High Compression
High-density 48,06 m³ / 18,14 KGF compressive