

## Serie V-1100

Firm mattress consisting of a lamination of high-resilience foam and soft foam ensuring a soft, uniform feeling of comfort. Supports a weight up to 700 lbs (317.5 kg). Easy maintenance, the Series V-1100 allows a variety of formats. It is treated antibacterial, antimold, antifungal and flame retardant. Flexible polyurethane-polyester fabric helps to reduce pressure.

The V-1100 Series is covered with an Integrity 30 fire net to make it self-extinguishing. It meets flammability standards CAN/ULC-S137:2017, CAL TB-129 and CAL 117.



Zippered cover made of antibacterial, anti-fungal, anti-mold and fire resistant polyurethane laminated polyester.

- 6" height
  - 5-year warranty  
(1 year 100% / pro-rated over 5 years)
- Shipping and handling cost are not included.



VISCO-ELASTIC FOAM  
High-density 4.0 lb ft<sup>3</sup> /  
11 ILD compressive

HIGH DENSITY HR FOAM  
ULTRA COMPRESSED HIGH COMPRESSION  
High-density 3.0 lb ft<sup>3</sup> / 40 lb ILD compressive