



MultiPositions



Bed esthetically designed
for home care



Made in Canada
Authentic canadian product.



Ultra-quiet DC motors
The quietest on the market,
keeps a calm environment.



Battery
for ensuring temporary autonomous.

Technical specifications

	36" (91 cm)
	39" (99 cm)
Platform width:	54" (137 cm)
	60" (152 cm)
Platform length:	80" (203 cm)
Minimum height:	14"/35.5 cm
Maximum height:	24"/60,9 cm
Load capacity:	500 lb/227 kg
Maximum backrest angle:	68.5° (degrees)
Maximum footrest angle:	38° (degrees)
Trendelenburg Position:	-10° à 10° (degrees)
	36"/91 cm = 235 lb/107 kg
	39"/99 cm = 260 lb/118 kg
Bed weight:	54"/137 cm = 295 lb/134 kg
	60"/152 cm = 310 lb/141 kg

Features

- Bed esthetically designed for home care
- Adjustable bed at the head and feet
- Height adjustable bed
- Ultra quiet DC motors. Helps maintain a calm and relaxed environment.
- Cordless remote control, easy to use
- Battery ensuring temporary autonomous operation
- Castors with brake
- Melamine bedspring 1"/2.54 cm
- Mattress retainer at the foot of the bed
- Authentic Canadian product
- CSA certified to Canadian standards

Options

- **Remote control:**
BT-COAR wired, digital clock and 2 memory position
RF-COAR wireless, digital clock and 2 memory position

Multipositions

Accessories:

- Decorative melamine headboard and footboard (color choice: white, rosewood, maple, and other colors on request)
- Rotec rotating half side rails
- Rotec rotating half side rails with increased protection height
- Rotec tilt-up half side rails 28" / 71.12 cm
- Rotec tilt-up half side rails 23" / 58.42 cm
- Metallic half side rails
- Side mattress retainer
- Foot mattress retainer
- Metallic support with bar and trapeze
- Stabilizer plate
- Clone bed
- Fix legs
- Side rail pad