

Series R-100

Medium-firm mattress consisting of a high-resiliency foam lamination and soft foam ensuring smooth, even comfort. The Series R-100 Model is low-maintenance and available in a variety of sizes. It undergoes antibacterial, anti-mold, anti-fungal and fire resistance treatment. The flexible polyurethane-polyester fabric helps reduce pressure.

The series R-100 is covered with an Integrity 30 fire barrier fabric to make it self-extinguishing.

Meet: CAN/ULC-S137:2017 standard
CAL TB 129 standard.




Zippered cover made of antibacterial, anti-fungal, anti-mold and fire resistant polyurethane laminated polyester.

Ultra-strong handles

- 6" height
- 5-year warranty
(1 year 100% / pro-rated over 5 years)

Shipping and handling cost are not included.



High-density 2.0 lb ft³ / 30 ILD compressive strength foam layer

High-density 2.0 lb ft³ / 55 lb ILD compressive strength foam side supports

High-density 1.5 lb ft³ / 40 lb ILD compressive strength foam layer

Standard sizes available:

35" x 80" x 6"	89 cm x 203 cm x 15.24 cm
36" x 80" x 6"	91.4 cm x 203 cm x 15.24 cm
39" x 80" x 6"	99 cm x 203 cm x 15.24 cm
48" x 80" x 6"	122 cm x 203 cm x 15.24 cm

Custom sizes (on request)