

Brilliantly designed

Series R-300

Soft mattress consisting of a high-resiliency foam center covered with memory foam to prevent pressure sores while providing smooth, even comfort. The Series R-300 model is low-maintenance and available in a variety of sizes. It is antibacterial, anti-mold, anti-fungal and fire resistant. The flexible polyurethane-polyester fabric helps reduce pressure.

The Series R-300 is covered with an Integrity 30 fire barrier fabric to make it self-extinguishing. Meet: CAN/ULC-S137:2017 standard
CAL TB 129 standard.




Zippered cover made of antibacterial, anti-fungal, anti-mold and fire resistant polyurethane laminated polyester.

Ultra-strong handles

- 6" height
- 5-year warranty
(1 year 100% / pro-rated over 5 years)

Shipping and handling cost are not included.



Serene™ quick release memory foam comfort layer
High-density 2.5 lb ft³ / 12 lb ILD compressive strength

High-density 2.0 lb ft³ / 70 lb ILD compressive strength foam side supports

High-density 2.5 lb ft³ / 40 lb ILD compressive strength foam core with vented channels

Standard sizes available:

35" x 80" x 6"	89 cm x 203 cm x 15.24 cm
36" x 80" x 6"	91.4 cm x 203 cm x 15.24 cm
39" x 80" x 6"	99 cm x 203 cm x 15.24 cm
48" x 80" x 6"	122 cm x 203 cm x 15.24 cm

Custom sizes (on request)

* Photos and specifications may change without prior notice.